
TECHNOPLASTIN® HIS

REF	5003009	TECHNOPLASTIN® HIS	12 x 2 mL
REF	5003010	TECHNOPLASTIN® HIS f. Ceveron	12 x 4 mL
REF	5003030	TECHNOPLASTIN® HIS	2 x 10 mL
REF	5003026	TECHNOPLASTIN® HIS	6 x 10 mL
REF	5003028	TECHNOPLASTIN® HIS f. Ceveron	6 x 20 mL
REF	5003021	TECHNOPLASTIN® HIS	20 x 10 mL
REF	5003022	TECHNOPLASTIN® HIS f. Ceveron	20 x 20 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Kľúčová slova / Značenje simbola			
	manufacturer / Hersteller / fabbricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производитель / Производитель / výrobc / Proizvođač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace / срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba fídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
	CE-mark / CE-Kennzeichnung / marchio CE / marca de CE / Simbolo da CE / marquage CE / CE-märkning / CE-mærket / CE-merke / CE-σημάδι / CE марка / CE-označení / маркировка CE / značka CE / CE-marka		determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelse / Bestemmelse / προσδιορισμοί / брой тестове / stanovení / определний / počet stanovení / Definicija
AQUA	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / деостилирана вода / destilovaná voda / дистилювана вода / destilovaná voda / Destilisaná Voda	LOT	lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / партиа / партида номер / šarže / лот / šarže / Serija
BUF	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionspuffer / Reaktionsbuffer / Reaktionsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer	MTP	microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροπιλοδότησης / Микротительна пласка / Микропланшет / Mikrotitracijska destička / Mikrotitrationsplatte
CAL	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	REF	catalogue number / Katalognummer / numero de catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
CONJ	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	RTU	ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / färdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
CONT	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контроль / Контрольный образец / Kontrola / Kontrola	STOP	stop solution / Stopplösung / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп раствор / Стоп-раствор / Zastavovací roztok / Stop solucija
DIL	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / späd eller upplös i / fortyndes eller opløses i / Fortyndes eller oppløses i / αραιωση ή διάλυση σε / разтворете или разрежете с / zředí nebo rozpustí v / разбавить или растворить в / naředí nebo rozpustí v / razrediti ili rastvoriti u	SUB	substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
INC	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationspuffer / Inkubationsbuffer / Inkubationsbuffer / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer	WASH	washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösingskoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миелс разтвор / Концентрат промывочного раствора / Koncentrat promývacieho roztoku / Koncentrat solucije za ispiranje
IVD	in vitro diagnostic use / in vitro Diagnostikum / diagnostico in vitro / diagnóstico in vitro / diagnóstico in vitro / diagnostico in vitro / diagnostico in vitro / diagnostico in vitro / διάλυμα για in vitro / для in vitro / для in vitro / diagnostický prostředek in vitro / in vitro diagnostika		

TECHNOPLASTIN® HIS

Ca-thromboplastin for the determination of the prothrombin time (PT, Quick)

GB

PRODUCT DESCRIPTION INTENDEND USE

TECHNOPLASTIN® HIS is used for the determination of the prothrombin time (Quick test). This screening test is used

- for the control of blood coagulation disorders of the extrinsic system,
- for the surveillance of oral anticoagulant therapy,
- for separate determination of individual coagulation factors of the extrinsic system,
- to check the liver's function of synthesizing coagulation factors in hepatopathies.

A heparin neutralizing agent in the reagent permits the assessment of plasma obtained from patients on conventional heparin therapy (0.2 - 0.8 U. hep./mL).

COMPOSITION

TECHNOPLASTIN® HIS (HIS = HeparinInSensitive) is a standardised Ca-thromboplastin reagent obtained from rabbit brain which is characterized by sensitivity to coagulation factors II, V, VII and X. In addition, the reagent contains a heparin neutralizing agent.

MATERIAL REQUIRED (not supplied with the kit)

- Pipettes
- Distilled water
- Control Plasmas and Calibrators

REF	5020040	Coagulation Control N **	5 x 1 mL
REF	5011050	Coagulation Control AK **	5 x 1 mL
REF	5010004	AK-Calibrant	4 x 1 mL

** or any other package sizes, special Ceveron® alpha or Technoclot Control and Calibration reagents of Technoclone.

WARNING AND PRECAUTIONS

- IVD for *in vitro* diagnostic use
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...+8 °C.

Stability after reconstitution:

37 °C 6 hours	RT* 24 hours	12 °C (Ceveron) 5 days	+2...+8 °C 6 days	-20 °C 2 months
------------------	-----------------	---------------------------	----------------------	--------------------

Upon storage, caps should be screwed tightly. * = room temperature

TEST PROCEDURE

PREPARATION OF PLASMA SAMPLES

Plasma separation: Mix 9 parts of venous blood and 1 part of sodium citrate solution (0.11 mol/L) and centrifuge for 15 min. at an RCF of at least 2500 (corresponding to DIN 58905). Keep the plasma at room temperature (up to 4 hours).

PREPARATION OF REAGENT

The lyophilized reagent is to be reconstituted with the indicated volume of distilled water at room temperature and may be used immediately. For Standardisation a reconstitution time of 30 min is recommendend.

PERFORMANCE OF THE TEST

CEVERON

Technoclone provides Application sheets for Ceveron® alpha. The Application sheets contain analyser/assay specific handling and performance information which may differ from that provided in these instruction for use. In this case the information contained in the Application sheets supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

MANUAL

Pipette the following substances into a test tube preheated to 37°C:

	Citrated plasma	0.1 mL
	Incubate for 1 min at + 37°C	
+	TECHNOPLASTIN® HIS (37°C)	0.2 mL

Upon addition of TECHNOPLASTIN® HIS, start the stopwatch or the coagulometer resp. and determine the clotting time.

ANALYSES RESULTS

CALCULATION OF THE RESULTS

The prothrombin time is indicated in seconds, in % of normal (Quick), or in INR. Values can be converted by means of the enclosed reference curve in form of a table, considering the indicated method. Reference curves have to be checked for accuracy by means of control plasma prior to use. If the values are beyond the confidence range given for control plasma, a separate reference curve has to be prepared.

For the preparation of the reference curve Technoclone's calibration plasma or the following dilutions of a normal plasma pool (Coagulation Control N) with isotonic saline solution may be used.

- % of normal dilution	100 % undil.	50 % 1+1	25 % 1+3	12.5 % 1+7	10 % 1+9
---------------------------	-----------------	-------------	-------------	---------------	-------------

These dilutions are to be determined like patients samples. The mean values of triple determinations are plotted on a reciprocal curve sheet and joined linearly.

The values indicated in % of normal can be transformed to INR according to the enclosed reference curve table. The indication of the prothrombin time in INR, as recommended by the WHO, is based on the following formula:

$$INR = \left(\frac{PT \text{ patient plasma (sec)}}{PT \text{ normal plasma (sec)}} \right)^{ISI}$$

REFERENCE RANGE

Normal time – approximate value: 13 – 17 sec.

Normal range in % of normal: 70 - > 100 %

Therapeutic range of oral anticoagulants: INR of 2.0 – 4.5 of 20 – 42% of normal, resp.

STANDARDISATION

The reference batch is calibrated against the Secondary Int. Reference Preparation for rabbit brain thromboplastin reagent (plain), RBT/79.

The standardisation is performed by use of the "Simplified calibration procedure" (Loeliger et al. Thrombosis Diathesis Haemorrh 33-172, 1975 and Loeliger et al. Thrombosis and Haemostasis 40,272, 1978).

LIMITATION OF THE TEST

Inhibitors of the Lupus type anticoagulant can influence prothrombin time and lead to INRs that do not accurately reflect the true level of anticoagulation.

The choice of anticoagulant (i.e. oxalate instead of citrate) and the condition of the specimen (e.g. hemolyzed, lipemic, parenteral feeding, etc.) may affect results. The latter is particularly true of optical instrumentation measurement of the PT.

PERFORMANCE CHARACTERISTICS

Performance data are given below. Results obtained in individual laboratories may differ.

PRECISION

Reproducibility was determined with different samples (in series and day to day). The following results were obtained:

Sample	Intra assay		Inter assay	
	Sample 1	Sample 2	Sample 1	Sample 2
n	12	12	24	24
MV INR	0,99	3,14	1,00	3,15
SD (%)	0,012	0,037	0,010	0,039
CV (%)	1,24	1,18	0,91	1,24

COMPARISON OF METHODS OR CORRELATION

Following correlation (INR) was obtained in comparing TECHNOPLASTIN HIS with :

Recombiplastin 2G	y = 1.073x	R ² = 0.9837
Thromborel S	y = 1.0082x	R ² = 0.9927

LITERATURE

Moll, St and Ortel, Th. L. Monitoring Warfarin Therapy in Patients with Lupus Anticoagulants. Ann Int Med 127, 177-185 (1997).

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