

Performance evaluation of automated C1 esterase measurement on the new coagulation analyzer Ceveron® c100

M. Unterberger, L. Wagner, N. B. Binder

Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH, Vienna, Austria

INTRODUCTION

Background: The most important laboratory parameter for correct diagnosis of hereditary angioedema (HAE) or angioedema due to acquired C1-Inhibitor (C1-INH) deficiency is reduced C1-INH function. Fully automated and precise measurement of C1-INH function is a valuable diagnostic cool.


Aim: Aim of the study was to evaluate assay performance of the chromogenic C1-INH assay on the new coagulation analyzer Ceveron[®] c100 in regards of accuracy and precision between runs.

MATERIALS AND METHODS

Method: To evaluate assay performance, lyophilized plasma samples were measured automated on the coagulation analyzer Ceveron[®] c100 with the chromogenic C1-INH assay kit. The assay is based on the inhibition of an excess of added Esterase by the C1-INH in the plasma sample, as measured by a chromogenic C1-INH substrate. The precision, accuracy, detection limit (LoD) and linearity of automated measurement of C1-INH were verified.



RESULTS



 $R^2 = 1.0 \pm 0.1$

Fig. 1 Calibration curve on Ceveron c100 analyzer was made in the range of 0-125 IU/dL using the standard analyzer settings for C1-INH calibration with a calibration plasma traceable to the WHO 1st International Standard for C1-Inhibitor, plasma.

Mean IU/dL C1-INH dean IU/dL C1-INH timepoint 2 % deviation

Coagulation Control N 1P72BC1K 128.8 125.9 -2.30

Coagulation Control A 3P72BC1K 54.0 51.6 -4.40

Tab. 1 To proof the stability of calibration curves, recovery of controls was calculated at different time points. Recovery of controls

The detection limit for C1-INH was determined with 6.60 IU/dL and lower limit of blanc with 2.2 IU/dL

Performance of C1-INH measurements

Ceveron® c100 Intra-assay CVs		Sample			
		2	3	4	
		n=11	n=11	n=11	
mean	89.0	53.6	78.7	71.3	
SD	4.1	3.8	4.4	4.8	
CV	4.57%	7.07%	5.62%	6.72%	
Ceveron® c100 Inter-assay CVs		nple			
		2 n=3			
	110.0	40 N			
mean	110.8	46.0			
	mean SD CV n® c100 say CVs	say CVs	1 2 n=11 n=11 mean 89.0 53.6 SD 4.1 3.8 CV 4.57% 7.07% Sample say CVs 1 2 n=3 n=3	1 2 3 n=11 n=11 n=11 mean 89.0 53.6 78.7 SD 4.1 3.8 4.4 CV 4.57% 7.07% 5.62% Sample say CVs 1 2 n=3 n=3	

Tab.	2a a	n 2b	Precision	was ve	ry good with
intra-	assay	and in	nter-assay	variatio	ons < 10%.

	140,0									
	120,0 -			-						
	5 100,0 -									
	Pa 80,0									
	Weasured IU/dl 80,0 - 60,0 - 40,0 -	مرمد مرمد مرمد مرمد مرمد مرمد مرمد مرمد								
	2 40,0 -	0								
	20,0									
	0,0									
	0 20 40 60 80 100 120 140 Calculated IU/dL									
	Linear r	egression	correlation							
	1/slope	R square	Correla- tion coeffi-	P value						
	0.9908	0.9944	0.993	< 0.0001						

Fig. 3 Dilution linearity was tested from 125 IU/dL C1-INH to 0.




Fig. 4 Recovery of controls was within 100% ± 10% of target value

CONCLUSIONS

Our data demonstrate that using the chromogenic assay kit TECHNOCHROM® C1-INH in optimized settings on Ceveron® c100 the determination of functional C1-INH can be performed with very good performance within less than 10 minutes.