

ISI-Werte / ISI - values

System	ISI
Ceveron one, four, ten	1.02
Thrombotrack - Coagulometer	1.04
Thrombotimer	1.05
Thrombolyzer	0.98

System	ISI	MNPT	sec	%
Ceveron m1, m2, m4	1.08	16.2	15.1	111
			30.6	44
			41.4	29.6
			56.6	21.3

TPZ Tabelle / PT Table: ● Ceveron one, four, ten

Pat.sec	13.0 sec.		13.5 sec.		14.0 sec.		14.5 sec.		15.0 sec.		15.5 sec.	
	INR	%	INR	%	INR	%	INR	%	INR	%	INR	%
11.0	-	133	-	142	-	153	-	164	-	177	-	191
11.5	-	122	-	131	-	140	-	150	-	161	-	172
12.0	-	114	-	121	-	129	-	138	-	147	-	157
12.5	-	106	-	113	-	120	-	128	-	136	-	145
13.0	1.00	100	-	106	-	113	-	120	-	127	-	135
13.5	1.04	94	1.00	100	-	106	-	112	-	119	-	126
14.0	1.08	89	1.04	95	1.00	100	-	106	-	112	-	118
14.5	1.12	85	1.08	90	1.04	95	1.00	100	-	106	-	111
15.0	1.16	81	1.11	85	1.07	90	1.04	95	1.00	100	-	105
15.5	1.20	77	1.15	82	1.11	86	1.07	90	1.03	95	1.00	100
16.0	1.23	74	1.19	78	1.15	82	1.11	86	1.07	91	1.03	95
16.5	1.27	71	1.23	75	1.18	79	1.14	83	1.10	87	1.07	91
17.0	1.31	68	1.26	72	1.22	76	1.18	79	1.14	83	1.10	87
18.0	1.39	64	1.34	67	1.29	70	1.25	73	1.20	77	1.16	80
19.0	1.47	59	1.42	62	1.36	65	1.32	68	1.27	71	1.23	75
20.0	1.55	56	1.49	58	1.44	61	1.39	64	1.34	67	1.30	70
21.0	1.63	53	1.57	55	1.51	58	1.46	60	1.41	63	1.36	65
22.0	1.71	50	1.64	52	1.58	54	1.53	57	1.48	59	1.43	62
23.0	1.79	47	1.72	50	1.66	52	1.60	54	1.54	56	1.49	58
24.0	1.86	45	1.79	47	1.73	49	1.67	51	1.61	53	1.56	55
25.0	1.94	43	1.87	45	1.80	47	1.74	49	1.68	51	1.63	53
26.0	2.02	41	1.95	43	1.88	45	1.81	47	1.75	49	1.69	50
27.0	2.10	40	2.02	41	1.95	43	1.88	45	1.82	47	1.76	48
28.0	2.18	38	2.10	40	2.02	41	1.95	43	1.89	45	1.82	46
29.0	2.26	37	2.17	38	2.10	40	2.02	41	1.95	43	1.89	45
31.0	2.42	34	2.33	36	2.24	37	2.16	38	2.09	40	2.02	41
33.0	2.58	32	2.48	33	2.39	35	2.31	36	2.23	37	2.16	39
35.0	2.74	30	2.63	31	2.54	33	2.45	34	2.37	35	2.29	36
37.0	2.89	29	2.79	30	2.68	31	2.59	32	2.50	33	2.42	34
39.0	3.05	27	2.94	28	2.83	29	2.73	30	2.64	31	2.55	32
41.0	3.21	26	3.09	27	2.98	28	2.88	29	2.78	30	2.69	31
43.0	3.37	25	3.25	26	3.13	26	3.02	27	2.92	28	2.82	29
45.0	3.53	24	3.40	24	3.28	25	3.16	26	3.05	27	2.95	28
47.0	3.69	23	3.55	23	3.42	24	3.30	25	3.19	26	3.09	27
49.0	3.85	22	3.71	22	3.57	23	3.45	24	3.33	25	3.22	26
51.0	4.01	21	3.86	22	3.72	22	3.59	23	3.47	24	3.35	25
53.0	4.17	20	4.01	21	3.87	22	3.73	22	3.61	23	3.49	24
55.0	4.33	19	4.17	20	4.02	21	3.88	22	3.74	22	3.62	23
59.0	4.65	18	4.48	19	4.31	19	4.16	20	4.02	21	3.89	21
63.0	4.97	17	4.78	18	4.61	18	4.45	19	4.30	20	4.16	20
67.0	5.29	16	5.09	17	4.91	17	4.74	18	4.58	18	4.43	19
71.0	5.61	15	5.40	16	5.21	16	5.02	17	4.85	17	4.69	18
75.0	5.93	15	5.71	15	5.50	16	5.31	16	5.13	17	4.96	17
80.0	6.34	14	6.10	14	5.88	15	5.67	15	5.48	16	5.30	16
85.0	6.74	13	6.49	13	6.25	14	6.03	14	5.83	15	5.64	15

TPZ Tabelle / PT Table: ● Thrombotrack - Coagulometer

Pat. sec	13.0 sec.		13.5 sec.		14.0 sec.		14.5 sec.		15.0 sec.		15.5 sec.	
	INR	%	INR	%	INR	%	INR	%	INR	%	INR	%
11.0	-	134	-	144	-	155	-	167	-	180	-	195
11.5	-	123	-	132	-	141	-	152	-	163	-	175
12.0	-	114	-	122	-	130	-	139	-	149	-	160
12.5	-	107	-	114	-	121	-	129	-	137	-	147
13.0	1.00	100	-	106	-	113	-	120	-	128	-	136
13.5	1.04	94	1.00	100	-	106	-	112	-	119	-	127
14.0	1.08	89	1.04	94	1.00	100	-	106	-	112	-	119
14.5	1.12	85	1.08	90	1.04	95	1.00	100	-	106	-	112
15.0	1.16	81	1.12	85	1.07	90	1.04	95	1.00	100	-	105
15.5	1.20	77	1.15	81	1.11	86	1.07	90	1.03	95	1.00	100
16.0	1.24	74	1.19	78	1.15	82	1.11	86	1.07	90	1.03	95
16.5	1.28	71	1.23	74	1.19	78	1.14	82	1.10	86	1.07	91
17.0	1.32	68	1.27	71	1.22	75	1.18	79	1.14	83	1.10	87
18.0	1.40	63	1.35	66	1.30	69	1.25	73	1.21	76	1.17	80
19.0	1.48	59	1.43	62	1.37	65	1.32	68	1.28	71	1.24	74
20.0	1.57	55	1.51	58	1.45	61	1.40	63	1.35	66	1.30	69
21.0	1.65	52	1.58	54	1.53	57	1.47	59	1.42	62	1.37	65
22.0	1.73	49	1.66	51	1.60	54	1.54	56	1.49	59	1.44	61
23.0	1.81	47	1.74	49	1.68	51	1.62	53	1.56	55	1.51	58
24.0	1.89	44	1.82	46	1.75	48	1.69	51	1.63	53	1.58	55
25.0	1.98	42	1.90	44	1.83	46	1.76	48	1.70	50	1.64	52
26.0	2.06	41	1.98	42	1.90	44	1.84	46	1.77	48	1.71	50
27.0	2.14	39	2.06	41	1.98	42	1.91	44	1.84	46	1.78	48
28.0	2.22	37	2.14	39	2.06	41	1.98	42	1.91	44	1.85	46
29.0	2.31	36	2.22	38	2.13	39	2.06	41	1.99	42	1.92	44
31.0	2.47	34	2.38	35	2.29	36	2.21	38	2.13	39	2.06	41
33.0	2.64	31	2.54	33	2.44	34	2.35	35	2.27	37	2.20	38
35.0	2.80	29	2.70	31	2.60	32	2.50	33	2.42	34	2.33	36
37.0	2.97	28	2.86	29	2.75	30	2.65	31	2.56	32	2.47	33
39.0	3.14	26	3.02	27	2.90	28	2.80	30	2.70	31	2.61	32
41.0	3.31	25	3.18	26	3.06	27	2.95	28	2.85	29	2.75	30
43.0	3.47	24	3.34	25	3.22	26	3.10	27	2.99	28	2.89	29
45.0	3.64	23	3.50	24	3.37	25	3.25	25	3.14	26	3.03	27
47.0	3.81	22	3.66	23	3.53	24	3.40	24	3.28	25	3.17	26
49.0	3.98	21	3.83	22	3.68	23	3.55	23	3.43	24	3.31	25
51.0	4.15	20	3.99	21	3.84	22	3.70	22	3.57	23	3.45	24
53.0	4.32	19	4.15	20	4.00	21	3.85	22	3.72	22	3.60	23
55.0	4.49	19	4.31	19	4.15	20	4.01	21	3.87	22	3.74	22
59.0	4.83	18	4.64	18	4.47	19	4.31	19	4.16	20	4.02	21
63.0	5.17	16	4.97	17	4.79	18	4.61	18	4.45	19	4.30	20
67.0	5.51	16	5.30	16	5.10	17	4.92	17	4.75	18	4.59	18
71.0	5.85	15	5.63	15	5.42	16	5.22	16	5.04	17	4.87	17
75.0	6.20	14	5.96	14	5.74	15	5.53	15	5.34	16	5.16	16
80.0	6.63	13	6.37	14	6.14	14	5.92	15	5.71	15	5.52	15
85.0	7.06	12	6.79	13	6.54	13	6.30	14	6.08	14	5.88	15

Pat. sec	13.0 sec.		13.5 sec.		14.0 sec.		14.5 sec.		15.0 sec.		15.5 sec.	
	INR	%	INR	%	INR	%	INR	%	INR	%	INR	%
11.0	-	134	-	144	-	155	-	167	-	181	-	195
11.5	-	123	-	132	-	142	-	152	-	163	-	176
12.0	-	114	-	122	-	130	-	140	-	149	-	160
12.5	-	107	-	114	-	121	-	129	-	138	-	147
13.0	1.00	100	-	106	-	113	-	120	-	128	-	136
13.5	1.04	94	1.00	100	-	106	-	113	-	119	-	127
14.0	1.08	89	1.04	94	1.00	100	-	106	-	112	-	119
14.5	1.12	85	1.08	90	1.04	95	1.00	100	-	106	-	112
15.0	1.16	81	1.12	85	1.07	90	1.04	95	1.00	100	-	105
15.5	1.20	77	1.16	81	1.11	86	1.07	90	1.03	95	1.00	100
16.0	1.24	74	1.19	78	1.15	82	1.11	86	1.07	90	1.03	95
16.5	1.28	71	1.23	74	1.19	78	1.14	82	1.10	86	1.07	91
17.0	1.32	68	1.27	71	1.23	75	1.18	79	1.14	83	1.10	87
18.0	1.41	63	1.35	66	1.30	69	1.25	73	1.21	76	1.17	80
19.0	1.49	59	1.43	62	1.38	65	1.33	68	1.28	71	1.24	74
20.0	1.57	55	1.51	58	1.45	60	1.40	63	1.35	66	1.31	69
21.0	1.65	52	1.59	54	1.53	57	1.47	59	1.42	62	1.37	65
22.0	1.73	49	1.67	51	1.60	54	1.55	56	1.49	58	1.44	61
23.0	1.82	47	1.75	49	1.68	51	1.62	53	1.56	55	1.51	58
24.0	1.90	44	1.82	46	1.76	48	1.69	50	1.63	52	1.58	55
25.0	1.98	42	1.90	44	1.83	46	1.77	48	1.71	50	1.65	52
26.0	2.06	40	1.98	42	1.91	44	1.84	46	1.78	48	1.72	50
27.0	2.15	39	2.06	40	1.99	42	1.92	44	1.85	46	1.79	47
28.0	2.23	37	2.14	39	2.06	40	1.99	42	1.92	44	1.86	45
29.0	2.31	36	2.22	37	2.14	39	2.06	40	1.99	42	1.92	44
31.0	2.48	33	2.38	35	2.30	36	2.21	38	2.14	39	2.06	40
33.0	2.65	31	2.55	33	2.45	34	2.36	35	2.28	36	2.20	38
35.0	2.82	29	2.71	31	2.61	32	2.51	33	2.42	34	2.34	35
37.0	2.98	28	2.87	29	2.76	30	2.66	31	2.57	32	2.48	33
39.0	3.15	26	3.03	27	2.92	28	2.81	29	2.71	30	2.62	32
41.0	3.32	25	3.19	26	3.07	27	2.96	28	2.86	29	2.76	30
43.0	3.49	24	3.36	25	3.23	26	3.12	27	3.01	28	2.91	28
45.0	3.66	23	3.52	24	3.39	24	3.27	25	3.15	26	3.05	27
47.0	3.83	22	3.68	23	3.55	23	3.42	24	3.30	25	3.19	26
49.0	4.00	21	3.85	22	3.70	22	3.57	23	3.45	24	3.33	25
51.0	4.17	20	4.01	21	3.86	22	3.72	22	3.59	23	3.47	24
53.0	4.34	19	4.18	20	4.02	21	3.88	22	3.74	22	3.61	23
55.0	4.52	19	4.34	19	4.18	20	4.03	21	3.89	21	3.76	22
59.0	4.86	17	4.67	18	4.50	19	4.34	19	4.18	20	4.04	21
63.0	5.20	16	5.00	17	4.82	18	4.64	18	4.48	19	4.33	19
67.0	5.55	15	5.34	16	5.14	17	4.95	17	4.78	18	4.62	18
71.0	5.90	15	5.67	15	5.46	16	5.26	16	5.08	17	4.91	17
75.0	6.25	14	6.00	14	5.78	15	5.57	15	5.38	16	5.20	16
80.0	6.68	13	6.42	14	6.18	14	5.96	14	5.75	15	5.56	15
85.0	7.12	12	6.84	13	6.59	13	6.35	14	6.13	14	5.92	15

Die Verwendbarkeit der Umrechnungstabellen muss durch eine Kontrolluntersuchung sichergestellt werden (z.B. mit Coagulation Control N von Technoclone). A control examination (e.g. by means of Technoclone's Coagulation Control N) must be performed to ensure the applicability of the conversion tables.

