



TECHNOPLASTIN HIS 10 mL

LOT 6T38CB0

3013942-161

ISI-Werte / ISI - values

System	ISI
Ceveron one, four, ten	1.08
Thrombotrack - Coagulometer	1.10
Thrombotimer	1.07
Thrombolyzer	1.02

System	ISI	MNPT	sec	%
Ceveron m1, m2, m4	1.06	15.5	14.4	111.0
			29.9	44.0
			40.0	29.6
			55.4	21.3

TPZ Tabelle / PT Table: ● Ceveron one, four, ten

Pat.sec	13.0 sec.		13.5 sec.		14.0 sec.		14.5 sec.		15.0 sec.		15.5 sec.	
	INR	%	INR	%	INR	%	INR	%	INR	%	INR	%
11.0	-	135	-	146	-	158	-	171	-	185	-	201
11.5	-	124	-	133	-	143	-	154	-	167	-	180
12.0	-	115	-	123	-	132	-	141	-	152	-	163
12.5	-	107	-	114	-	122	-	130	-	139	-	149
13.0	1.00	100	-	107	-	114	-	121	-	129	-	138
13.5	1.04	94	1.00	100	-	106	-	113	-	120	-	128
14.0	1.08	89	1.04	94	1.00	100	-	106	-	113	-	119
14.5	1.13	84	1.08	89	1.04	94	1.00	100	-	106	-	112
15.0	1.17	80	1.12	85	1.08	90	1.04	95	1.00	100	-	106
15.5	1.21	76	1.16	81	1.12	85	1.07	90	1.04	95	1.00	100
16.0	1.25	73	1.20	77	1.16	81	1.11	86	1.07	90	1.03	95
16.5	1.29	70	1.24	74	1.19	78	1.15	82	1.11	86	1.07	90
17.0	1.34	67	1.28	71	1.23	74	1.19	78	1.14	82	1.10	86
18.0	1.42	62	1.36	65	1.31	69	1.26	72	1.22	76	1.18	79
19.0	1.51	58	1.45	61	1.39	64	1.34	67	1.29	70	1.25	73
20.0	1.59	54	1.53	57	1.47	59	1.42	62	1.36	65	1.32	68
21.0	1.68	51	1.61	53	1.55	56	1.49	58	1.44	61	1.39	64
22.0	1.77	48	1.70	50	1.63	53	1.57	55	1.51	57	1.46	60
23.0	1.85	46	1.78	48	1.71	50	1.65	52	1.59	54	1.53	57
24.0	1.94	43	1.86	45	1.79	47	1.72	49	1.66	51	1.60	54
25.0	2.03	41	1.95	43	1.87	45	1.80	47	1.74	49	1.68	51
26.0	2.11	39	2.03	41	1.95	43	1.88	45	1.81	47	1.75	49
27.0	2.20	38	2.11	39	2.03	41	1.96	43	1.89	45	1.82	46
28.0	2.29	36	2.20	38	2.11	39	2.04	41	1.96	43	1.89	44
29.0	2.38	35	2.28	36	2.20	38	2.11	39	2.04	41	1.97	43
31.0	2.56	32	2.46	34	2.36	35	2.27	37	2.19	38	2.11	39
33.0	2.74	30	2.63	31	2.53	33	2.43	34	2.34	35	2.26	37
35.0	2.92	28	2.80	30	2.69	31	2.59	32	2.50	33	2.41	34
37.0	3.10	27	2.97	28	2.86	29	2.75	30	2.65	31	2.56	32
39.0	3.28	25	3.15	26	3.03	27	2.91	28	2.81	29	2.71	31
41.0	3.46	24	3.32	25	3.19	26	3.07	27	2.96	28	2.86	29
43.0	3.64	23	3.50	24	3.36	25	3.24	26	3.12	27	3.01	27
45.0	3.83	22	3.67	23	3.53	24	3.40	24	3.28	25	3.16	26
47.0	4.01	21	3.85	22	3.70	22	3.56	23	3.44	24	3.32	25
49.0	4.19	20	4.03	21	3.87	22	3.73	22	3.59	23	3.47	24
51.0	4.38	19	4.20	20	4.04	21	3.89	21	3.75	22	3.62	23
53.0	4.57	18	4.38	19	4.21	20	4.06	21	3.91	21	3.78	22
55.0	4.75	18	4.56	18	4.39	19	4.22	20	4.07	21	3.93	21
59.0	5.13	17	4.92	17	4.73	18	4.56	18	4.39	19	4.24	20
63.0	5.50	16	5.28	16	5.08	17	4.89	17	4.71	18	4.55	19
67.0	5.88	15	5.65	15	5.43	16	5.23	16	5.04	17	4.86	17
71.0	6.26	14	6.01	14	5.78	15	5.57	15	5.37	16	5.18	16
75.0	6.64	13	6.38	14	6.13	14	5.90	15	5.69	15	5.49	16
80.0	7.12	12	6.84	13	6.58	13	6.33	14	6.10	14	5.89	15
85.0	7.61	12	7.30	12	7.02	12	6.76	13	6.52	13	6.29	14

TPZ Tabelle / PT Table: ● Thrombotrack - Coagulometer

Pat. sec	13.0 sec.		13.5 sec.		14.0 sec.		14.5 sec.		15.0 sec.		15.5 sec.	
	INR	%	INR	%	INR	%	INR	%	INR	%	INR	%
11.0	-	136	-	147	-	159	-	173	-	188	-	205
11.5	-	125	-	134	-	145	-	156	-	169	-	183
12.0	-	115	-	124	-	133	-	142	-	153	-	165
12.5	-	107	-	114	-	122	-	131	-	140	-	151
13.0	1.00	100	-	107	-	114	-	122	-	130	-	139
13.5	1.04	94	1.00	100	-	106	-	113	-	121	-	128
14.0	1.09	89	1.04	94	1.00	100	-	106	-	113	-	120
14.5	1.13	84	1.08	89	1.04	94	1.00	100	-	106	-	112
15.0	1.17	80	1.12	84	1.08	89	1.04	95	1.00	100	-	106
15.5	1.21	76	1.16	80	1.12	85	1.08	90	1.04	95	1.00	100
16.0	1.26	72	1.21	77	1.16	81	1.11	85	1.07	90	1.04	95
16.5	1.30	69	1.25	73	1.20	77	1.15	81	1.11	86	1.07	90
17.0	1.34	66	1.29	70	1.24	74	1.19	78	1.15	82	1.11	86
18.0	1.43	61	1.37	65	1.32	68	1.27	72	1.22	75	1.18	79
19.0	1.52	57	1.46	60	1.40	63	1.35	66	1.30	69	1.25	73
20.0	1.61	53	1.54	56	1.48	59	1.43	62	1.37	65	1.32	68
21.0	1.70	50	1.63	53	1.56	55	1.50	58	1.45	61	1.40	63
22.0	1.79	47	1.71	50	1.65	52	1.58	54	1.53	57	1.47	59
23.0	1.88	45	1.80	47	1.73	49	1.66	51	1.60	54	1.55	56
24.0	1.97	43	1.89	45	1.81	47	1.74	49	1.68	51	1.62	53
25.0	2.06	41	1.97	42	1.90	44	1.82	46	1.76	48	1.69	50
26.0	2.15	39	2.06	41	1.98	42	1.90	44	1.83	46	1.77	48
27.0	2.24	37	2.15	39	2.06	40	1.99	42	1.91	44	1.84	46
28.0	2.33	36	2.24	37	2.15	39	2.07	40	1.99	42	1.92	44
29.0	2.42	34	2.33	36	2.23	37	2.15	39	2.07	40	2.00	42
31.0	2.61	32	2.50	33	2.40	34	2.31	36	2.23	37	2.15	39
33.0	2.80	30	2.68	31	2.58	32	2.48	33	2.39	35	2.30	36
35.0	2.98	28	2.86	29	2.75	30	2.64	31	2.55	32	2.46	34
37.0	3.17	26	3.04	27	2.92	28	2.81	29	2.71	31	2.61	32
39.0	3.36	25	3.22	26	3.10	27	2.98	28	2.87	29	2.77	30
41.0	3.55	23	3.41	24	3.27	25	3.15	26	3.03	27	2.93	28
43.0	3.74	22	3.59	23	3.45	24	3.32	25	3.20	26	3.08	27
45.0	3.94	21	3.78	22	3.63	23	3.49	24	3.36	25	3.24	26
47.0	4.13	20	3.96	21	3.81	22	3.66	23	3.53	24	3.40	24
49.0	4.32	19	4.15	20	3.98	21	3.83	22	3.69	23	3.56	23
51.0	4.52	19	4.34	19	4.16	20	4.01	21	3.86	22	3.72	22
53.0	4.72	18	4.52	19	4.35	19	4.18	20	4.03	21	3.88	21
55.0	4.91	17	4.71	18	4.53	19	4.35	19	4.19	20	4.05	21
59.0	5.31	16	5.09	17	4.89	17	4.71	18	4.53	19	4.37	19
63.0	5.71	15	5.47	16	5.26	16	5.06	17	4.87	17	4.70	18
67.0	6.11	14	5.86	15	5.63	15	5.41	16	5.22	16	5.03	17
71.0	6.51	13	6.25	14	6.00	14	5.77	15	5.56	15	5.36	16
75.0	6.92	13	6.64	13	6.37	14	6.13	14	5.91	15	5.70	15
80.0	7.43	12	7.12	12	6.84	13	6.58	13	6.34	14	6.12	14
85.0	7.94	11	7.62	12	7.32	12	7.04	12	6.78	13	6.54	13

TPZ Tabelle / PT Table: ● Thrombotimer

Pat. sec	13.0 sec.		13.5 sec.		14.0 sec.		14.5 sec.		15.0 sec.		15.5 sec.	
	INR	%	INR	%	INR	%	INR	%	INR	%	INR	%
11.0	-	135	-	145	-	157	-	170	-	184	-	199
11.5	-	124	-	133	-	143	-	154	-	166	-	179
12.0	-	115	-	123	-	131	-	141	-	151	-	162
12.5	-	107	-	114	-	122	-	130	-	139	-	148
13.0	1.00	100	-	107	-	113	-	121	-	129	-	137
13.5	1.04	94	1.00	100	-	106	-	113	-	120	-	127
14.0	1.08	89	1.04	94	1.00	100	-	106	-	112	-	119
14.5	1.12	84	1.08	89	1.04	95	1.00	100	-	106	-	112
15.0	1.17	80	1.12	85	1.08	90	1.04	95	1.00	100	-	106
15.5	1.21	76	1.16	81	1.12	85	1.07	90	1.04	95	1.00	100
16.0	1.25	73	1.20	77	1.15	81	1.11	86	1.07	90	1.03	95
16.5	1.29	70	1.24	74	1.19	78	1.15	82	1.11	86	1.07	91
17.0	1.33	67	1.28	71	1.23	74	1.19	78	1.14	82	1.10	86
18.0	1.42	62	1.36	65	1.31	69	1.26	72	1.22	76	1.17	79
19.0	1.50	58	1.44	61	1.39	64	1.34	67	1.29	70	1.24	73
20.0	1.59	54	1.52	57	1.46	60	1.41	63	1.36	65	1.31	68
21.0	1.67	51	1.60	54	1.54	56	1.49	59	1.43	61	1.38	64
22.0	1.76	48	1.69	51	1.62	53	1.56	55	1.51	58	1.45	60
23.0	1.84	46	1.77	48	1.70	50	1.64	52	1.58	55	1.53	57
24.0	1.93	44	1.85	46	1.78	48	1.71	50	1.65	52	1.60	54
25.0	2.01	42	1.93	43	1.86	45	1.79	47	1.73	49	1.67	51
26.0	2.10	40	2.02	41	1.94	43	1.87	45	1.80	47	1.74	49
27.0	2.19	38	2.10	40	2.02	41	1.95	43	1.88	45	1.81	47
28.0	2.27	37	2.18	38	2.10	40	2.02	41	1.95	43	1.88	45
29.0	2.36	35	2.27	37	2.18	38	2.10	40	2.02	41	1.96	43
31.0	2.53	33	2.43	34	2.34	35	2.26	37	2.17	38	2.10	40
33.0	2.71	31	2.60	32	2.50	33	2.41	34	2.33	36	2.25	37
35.0	2.89	29	2.77	30	2.67	31	2.57	32	2.48	33	2.39	35
37.0	3.06	27	2.94	28	2.83	29	2.73	30	2.63	31	2.54	33
39.0	3.24	26	3.11	27	2.99	28	2.88	29	2.78	30	2.68	31
41.0	3.42	24	3.28	25	3.16	26	3.04	27	2.93	28	2.83	29
43.0	3.60	23	3.46	24	3.32	25	3.20	26	3.09	27	2.98	28
45.0	3.78	22	3.63	23	3.49	24	3.36	25	3.24	26	3.13	26
47.0	3.96	21	3.80	22	3.66	23	3.52	24	3.40	24	3.28	25
49.0	4.14	20	3.97	21	3.82	22	3.68	23	3.55	23	3.43	24
51.0	4.32	19	4.15	20	3.99	21	3.84	22	3.71	22	3.58	23
53.0	4.50	19	4.32	19	4.16	20	4.00	21	3.86	22	3.73	22
55.0	4.68	18	4.50	19	4.32	19	4.17	20	4.02	21	3.88	22
59.0	5.05	17	4.85	17	4.66	18	4.49	19	4.33	19	4.18	20
63.0	5.41	16	5.20	16	5.00	17	4.82	18	4.65	18	4.49	19
67.0	5.78	15	5.55	15	5.34	16	5.15	17	4.96	17	4.79	18
71.0	6.15	14	5.91	15	5.68	15	5.47	16	5.28	16	5.10	17
75.0	6.53	13	6.27	14	6.03	14	5.81	15	5.60	15	5.41	16
80.0	6.99	13	6.71	13	6.46	13	6.22	14	6.00	14	5.79	15
85.0	7.46	12	7.17	12	6.89	13	6.64	13	6.40	14	6.18	14

Die Verwendbarkeit der Umrechnungstabellen muss durch eine Kontrolluntersuchung sichergestellt werden (z.B. mit Coagulation Control N von Technoclone). A control examination (e.g. by means of Technoclone's Coagulation Control N) must be performed to ensure the applicability of the conversion tables.

